PRODUCT SHEET

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- (1) Coma et al. "Automated separation of 99Tc using plastic scintillation resin PSresin and openview automated modular separation system (OPENVIEW-AMSS)", Journal of Radioanalytical and Nuclear Chemistry (2019) 321:1057–1065.
<https://doi.org/10.1007/s10967-019-06659-7>
- (2) Barrera et al. "A new plastic scintillation resin for single-step separation, concentration and measurement of technetium-99" Analytica Chimica Acta 936 (2016) 259-266.
<https://doi.org/10.1016/j.aca.2016.07.008>
- (3) Tarancón et al. "A new plastic scintillation resin for single-step separation, concentration and measurement of 99Tc", presented at the NRC9 (29/08/16 – 2/09/16, Helsinki, Finland)
- (4) Hidex eBook "Liquid Scintillation Measuring Procedures: New Developments"
<https://hidex.com/ebooks/liquid-scintillation-measuring-procedures/measuring-procedures/radionuclides-from-nuclear-fission-activities/2-3-14-tc-by-rad-disk-and-psresin/>
- (5) J. Garcia & A. Tarancon, "Radionuclide determinations with PS Resin MASS WaterRadd", presented at the European Users Group Meeting in Cambridge (UK) - 21/09/2018,
https://www.triskem-international.com/scripts/files/5bae2550c30ed4.50583030/11_j-garcia_a-tarancon_radionuclide-determinations-with-ps-resin_mass_waterradd.pdf
- (6) Giménez et al. "A new method based on selective fluorescent polymers (PSresin) for the analysis of 90Sr in presence of 210Pb in environmental samples" Applied Radiation and Isotopes, 199, 2023,
<https://doi.org/10.1016/j.apradiso.2023.110879>
- (7) Tarancón et al. "Recent applications of Plastic Scintillation Resins », presented at Raddec-Triskem Workshop 2024, (18/04/24, Portsmouth, UK)